




Listeria spp.
White colonies



L. monocytogenes
Blue colonies

The Difference Is Clear

Detect *Listeria monocytogenes* in 24 hours with the RAPID'L.Mono™ detection system: Reliable results that you can see.

Fast

- 24 hour pre-enrichment only
- Complete results in 48 hours
- Easy-to-read chromogenic reaction

Simple and Affordable

- Now available in bottled form
- No selective enrichment required
- Tedious biochemical tests are not necessary
- Rapid method at traditional price



Results

- Blue colonies are *L. monocytogenes*
- Green colonies with distinct yellow halo are *L. ivanovii*
- White colonies are other *Listeria* species

Validation

- AFNOR approved as confirmatory test
- Evaluation completed by USDA-ARS

For more information, visit www.bio-rad.com/ad/RapidLmono/




Visit us on the Web at discover.bio-rad.com
Call toll free at 1-800-4BIORAD (1-800-424-6723);
outside the US, contact your local sales office.



Journal of VETERINARY DIAGNOSTIC INVESTIGATION

Official Publication of the American Association of Veterinary Laboratory Diagnosticians, Inc.



Vol. 18 No. 2

Journal of Veterinary Diagnostic Investigation

March 2006 Pages 139-232