

# The Difference Is Clear

**Detect *Listeria monocytogenes* in 24 hours with the RAPID'L.Mono™ detection system: Reliable results that you can see.**

### Fast

- 24 hour pre-enrichment only
- Complete results in 48 hours
- Easy-to-read chromogenic reaction

### Simple and Affordable

- Now available in bottled form
- No selective enrichment required
- Tedious biochemical tests are not necessary
- Rapid method at traditional price


### Results

- Blue colonies are *L. monocytogenes*
- Green colonies with distinct yellow halo are *L. ivanovii*
- White colonies are other *Listeria* species

### Validation

- AFNOR approved as confirmatory test
- Evaluation completed by USDA-ARS

For more information, visit [www.bio-rad.com/ad/RapidLmono/](http://www.bio-rad.com/ad/RapidLmono/)




Strains of RAPID'L.Mono were independently evaluated by the AOAC Research Institute and were found to perform to the producer's specifications as listed in the test kit's description sheet. The producer certifies its conformity to all aspects of the AOAC International Performance Tested Certificate number 000400.

Visit us on the Web at [discover.bio-rad.com](http://discover.bio-rad.com)  
Call toll free at 1-800-4BIORAD (1-800-424-6723);  
outside the US, contact your local sales office.



# Journal of VETERINARY DIAGNOSTIC INVESTIGATION

Official Publication of the American Association of Veterinary Laboratory Diagnosticians, Inc.



Vol. 18 No. 3

Journal of Veterinary Diagnostic Investigation

May 2006 Pages 233-318